



NHS

Blood and Transplant

How to use RCI Assist Referral Support Tool

Helen Thom

RCI Development Lead – Transfusion 2024

Case study examples

- Please use the slides to learn how to navigate RCI Assist
- 5 case examples with panels have been provided
- Work your way through each case to practice using the tool
- The outcome will vary depending on what reagents you have available in your laboratory
- If you require any assistance or advice, please contact helen.thom@nhsbt.nhs.uk

RCI Assist



Blood and Transplant

What is RCI Assist?

As part of the integrated transfusion services strategic plan for Transfusion 2024, a referral support tool, RCI Assist, has been developed.


RCI Assist guides users when processing samples in the Hospital Transfusion Laboratory (HTL). RCI Assist contains decision points to either refer cases to Red Cell Immunohaematology (RCI) or resolve in-house and is supported with pop up information to improve HTL staff confidence.

Pathways include:

- Grouping
- Antibody Identification
- Antenatal
- Panreactives
- Elution
- Transfusion



Scan here to see how to navigate RCI Assist



How it works:




To move through the pathways, select the yes or no options next to the question boxes



To jump straight to a section, select the required tab



To get supporting information about a topic, click the information icon




Reach a decision point and either refer to RCI or result in-house

Accessing RCI Assist:



RCI Assist is available as a link in Sp-ICE and OBOS



Log in with your usual credentials



For further information on RCI Assist, please visit the [NBTC website](#) or email transfusion2024@nhsbt.nhs.uk

Example A

- 62-year-old female
- Group & Screen
- 2 units RBC
- “Anaemia”
- Not previously seen by your lab, but transfused at another hospital 5 years ago
- No historical antibodies
- Not previously seen by RCI



Patient group & screen results



Is ABO and D group interpretable?



YES

NO





Is ABO and D group interpretable?

NO



[Click to learn more](#)

Is ABO and D group interpretable?


Forward and reverse ABO groups match in accordance with Landsteiner's Law.

Note: Neonate groups are considered 'interpretable' from forward group only.



ABO forward group and D group show strong positive (3+/4+) reactions. Reverse group reactions may be weaker.



"For whichever ABO antigen is not present on the red cells, the corresponding antibody is found in the plasma"




Is the ABO and D group interpretable?



Is the ABO and D group interpretable?

The forward and reverse ABO group match in accordance with Landsteiner's law and the D group shows strong positive reactions






Is ABO and D group interpretable?



YES

NO





Is ABO and D group interpretable?



YES



Go to Antibody ID section

NO





Is antibody screen positive?



YES

NO



Is antibody screen positive?


NO





[Click to learn more](#)

Is antibody screen positive?

Positive result obtained by Indirect Antiglobulin Test (IAT) antibody screen.




Is the antibody screen positive?



Is the antibody screen positive?

Positive result obtained by IAT



	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other
1	R ₁ ^w R ₁	+	+	0	0	+	+	+	+	+	0	+	+	+	0	0	+	+	+	0	+	3
2	R ₂ R ₂	0	+	+	+	0	0	0	+	0	+	+	0	+	+	0	+	0	+	0	+	0
3	rr	0	0	0	+	+	0	+	0	0	+	0	0	+	0	+	0	+	0	+	0	3






Is antibody screen positive?



YES

NO





Is antibody screen positive?



YES



Check Patient History

Does patient have a known antibody?



YES

NO

Refer to RCI
(5)



NO





Is antibody screen positive?



YES

NO

Refer to RCI
(5)




Check Patient History

Does patient have a known
antibody?



YES

NO





Is antibody screen positive?



NO

YES

Refer to RCI
(5)



Check Patient History

Does patient have a known
antibody?




NO

[Click to learn more](#)

Check Patient History
Does patient have a known antibody?

Has this patient previously been tested and
found to have red cell antibodies?



62-year-old female

Group & Screen

2 units RBC

“Anaemia”

Not previously seen by your lab,
but transfused at another
hospital 5 years ago

No historical antibodies

Not previously seen by RCI





Is antibody screen positive?



YES

Check Patient History

Does patient have a known antibody?



YES

NO

Refer to RCI
(5)



NO




Help




Print page



Disclaimer: [Read Here](#)





[Click to learn more](#)

Perform antibody ID on primary panel


Perform the primary antibody identification panel used in your laboratory in accordance with manufactures' guidelines and your local procedures.





Antibody identification






Cell	Rh	D	C	E	c	e	M	N	S	s	Pl	Lu ^a	K	k	Kp^a	Le^a	Le^b	Fy^a	Fy^b	Jk^a	Jk ^b	Other	IAT	ENZ
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	0	0	0	+	0	0	+	+	0	+	0		2	0
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+	0	+	0	0	+	0	+		3	3
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	0	+	0	0	+	+	+		2	0
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+	0		0	0
5	r''r	0	0	+	+	+	0	+	0	+	4	0	0	+	0	0	+	+	0	+	+		0	0
6	rr	0	0	0	+	+	+	0	+	0	4	+	+	0	0	0	+	+	0	0	+		3	3
7	rr	0	0	0	+	+	0	+	0	+	2	0	+	+	0	+	0	0	+	+	0		3	3
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	+	0	+	0	+	+	+		0	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	0	+	0	+	0	+	0		2	0
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	0	+	0	+	0	+	0		0	0
<div>Probable anti-K and anti-M Have not excluded anti-Jkb and Lua</div>																						Auto	0	/
																						+ control	2	/





Possible anti-M and anti-K





Is antibody screen positive?

NO

YES

Check Patient History

Does patient have a known antibody?

YES

Refer to RCI
(5)

NO

Perform Antibody ID on
Primary Panel

Has a pan-reactive antibody
been identified by IAT?

NO

NO to 1 or

Check Patient History

1. Have antibody(ies) been identified?
2. Have additional clinically significant antibodies been excluded?
3. Is patient pregnant?
4. Has this previously been investigated by RCI?

YES to 1, 2 and 3

YES to 1 and

[Click to learn more](#)

Check Patient Results and History

1. Have antibody(ies) been identified?
2. Have additional clinically significant antibodies been excluded?
3. Is patient pregnant?
4. Has this previously been investigated by RCI?

Hints and tips – reading and interpreting antibody ID panels

- Know the clinical significance of antibodies (read BSH guidelines)
- Confirm antibodies as alloantibodies (antigen phenotyping)
- Don't forget antigen dosage of cells
- Remember which cells are destroyed by enzyme which is worth considering when you are doing exclusions
- Cold reactive alloantibodies may react by IAT (E.g. Lewis and M antibodies)
- Column Agglutination Technology (CAT) can be prone to non-specific reactions
- Sometimes enzyme pan-reactivity causes non-specific reactivity in IAT. In such cases, please ensure that you can exclude all clinically significant antibodies by IAT.

Antibody Inclusion



- Consist of red cells from eight or more group O donors
- For commonly encountered clinically significant alloantibodies: **2 antigen pos** and **2 antigen neg** cells
- 1 x R1R1 (CCDDee) and 1 x R1wR1 (CwCDDee)
- Between them, these two cells should express the antigens K, k, Fya, Fyb, Jka, Jkb, S, s
- One example of each of the phenotypes R2R2 (ccDDEE), r'r (Ccddee) and r''r (ccddEe)
- At least three examples of the phenotype rr (ccddee), including at least one K+, and collectively, homozygous expression of k, Jka, Jkb, S, s, Fya, and Fyb

Antibody Exclusion


- 6.2.5. When one antibody specificity has been identified, it is essential that the presence or absence of additional clinically significant antibodies is established
- 6.2.6. Failure to recognise all of the antibody specificities within a sample may lead to a haemolytic transfusion reaction (SHOT, [1996–2010](#)). In particular, the presence of **anti-Jk^a**, **anti-Jk^b**, **anti-S**, **anti-s**, **anti-Fy^a** and **anti-Fy^b** should be excluded using red cells having **homozygous** expression of the relevant antigen
- A single example only of each phenotype is sufficient for exclusion



Antibody identification



Cell	Rh	D	C	E	c	e	M	N	S	s	Pl	Lu ^a	K	k	Kp^a	Le^a	Le^b	Fy^a	Fy^b	Jk^a	Jk ^b	Other	IAT	ENZ
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	0	0	0	+	0	0	+	+	0	+	0		2	0
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+	0	+	0	0	+	0	+		3	3
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	0	+	0	0	+	+	+		2	0
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+	0		0	0
5	r''r	0	0	+	+	+	0	+	0	+	4	0	0	+	0	0	+	+	0	+	+		0	0
6	rr	0	0	0	+	+	+	0	+	0	4	+	+	0	0	0	+	+	0	0	+		3	3
7	rr	0	0	0	+	+	0	+	0	+	2	0	+	+	0	+	0	0	+	+	0		3	3
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	+	0	+	0	+	+	+		0	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	0	+	0	+	0	+	0		2	0
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	0	+	0	+	0	+	0		0	0
<div>Probable anti-K and anti-M Have not excluded anti-Jkb and Lua</div>																						Auto	0	/
																						+ control	2	/



What would you do next?



Good suggestions!



NHS

Blood and Transplant



Antibody screen

	Rh	C	D	E	c	e	C^w	M	N	S	s	P1	K	k	Kp^a	Le^a	Le^b	Fy^a	Fy^b	Jk^a	Jk^b	IAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	+	+	+	0	+	+	+	0	0	+	+	+	0	+	3
2	R ₂ R ₂	0	+	+	+	0	0	0	+	0	+	+	0	+	+	0	+	0	+	0	+	0
3	rr	0	0	0	+	+	0	+	0	0	+	0	0	+	0	+	0	+	0	+	0	3



Have now excluded anti-Jkb


What about anti-Lua?

Is anti-Lua clinically significant?



NHS

Blood and Transplant






BSH Guidelines – clinical significance of red cell antibodies

System	Specificity	Likely clinical significance in transfusion
Kell	Anti-K	Yes
MNS	Anti-M (active at 37°C)	Yes
Lu	Anti-Lu ^a	No

Antibody identification



Cell	Rh	D	C	E	c	e	M	N	S	s	Pi	Lu ^a	K	k	Kp ^a	L ^a	L ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT	ENZ
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	0	0	0	+	0	0	+	+	0	+	0	→	2	0
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+	0	+	0	0	+	0	+	→	3	3
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	0	+	0	0	+	+	+	→	2	0
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+	0	→	0	0
5	r''r	0	0	+	+	+	0	+	0	+	4	0	0	+	0	0	+	+	0	+	+	→	0	0
6	rr	0	0	→	→	→	+	0	+	→	+	+	0	0	0	+	+	0	0	0	+	→	3	3
7	rr	0	0	0	+	+	0	+	0	+	2	0	+	+	0	+	0	0	+	+	0	→	3	3
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	+	0	+	0	+	+	+	→	0	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	0	+	0	+	0	+	0		2	0
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	0	+	0	+	0	+	0	→	0	0
																						→	0	/
																						→	2	/



Anti-K and anti-M
Not excluded anti-Lua

Check Patient Results and History

- ✓ 1. Have antibody(ies) been identified?
- ✓ 2. Have additional clinically significant antibodies been excluded?
- ✗ 3. Is patient pregnant?
- ✗ 4. Has this previously been investigated by RCI?

Hints and tips – reading and interpreting antibody ID panels

- Know the clinical significance of antibodies (read BSH guidelines)
- Confirm antibodies as alloantibodies (antigen phenotyping)
- Don't forget antigen dosage of cells
- Remember which cells are destroyed by enzyme which is worth considering when you are doing exclusions
- Cold reactive alloantibodies may react by IAT (E.g. Lewis and M antibodies)
- Column Agglutination Technology (CAT) can be prone to non-specific reactions
- Sometimes enzyme pan-reactivity causes non-specific reactivity in IAT. In such cases, please ensure that you can exclude all clinically significant antibodies by IAT.

62-year-old female

Group & Screen

2 units RBC


“Anaemia”


Not previously seen by your lab,
but transfused at another
hospital 5 years ago


No historical antibodies

Not previously seen by RCI









62-year-old female

Group & Screen


2 units RBC


“Anaemia”

Not previously seen by your lab, but transfused at another hospital 5 years ago

No historical antibodies

Not previously seen by RCI







Is antibody screen positive?



NO

YES

Refer to RCI
(5)

Check Patient History

Does patient have a known
antibody?



YES

NO

Perform Antibody ID on
Primary Panel



Has a pan-reactive antibody
been identified by IAT?

YES

NO

NO to 1 or 2

Check Patient History



1. Have antibody(ies) been identified?
2. Have additional clinically significant antibodies been excluded?
3. Is patient pregnant?
4. Has this previously been investigated by RCI?

YES to 1, 2 and 3

YES to 1 and 2

Check Results



Does patient have Rh
system antibody with
positive auto? Consider
patient ethnicity and
IBGRL referral

Refer to RCI
(18i)

Optional
retrospective
referral (if YES to
1&2 but NO to 4)
Refer to RCI (8)



Has patient been
transfused in last
3 months?

YES

NO

Perform phenotype
for identified
antibody-antigen pair
if not already known



Refer to RCI (9)



Is patient pregnant
with an antibody
requiring quantification/
titration?



YES

NO



Is antibody screen positive?



NO

YES

Refer to RCI
(5)

Check Patient History

Does patient have
antibody[Click to learn more](#)

NO

Perform Antibody ID on
Primary PanelHas a pan-reactive antibody
been identified by IAT?

YES

NO

NO to 1 or 2

Check Patient History



1. Have antibody(ies) been identified?
2. Have additional clinically significant antibodies been excluded?
3. Is patient pregnant?
4. Has this previously been investigated by RCI?

YES to 1, 2 and 3

YES to 1 and 2

Check Results

Does patient have Rh
system antibody with
positive auto? Consider
patient ethnicity and
IBGRL referralRefer to RCI
(18i)Optional
retrospective
referral (if YES to
1&2 but NO to 4)
Refer to RCI (8)Has patient been
transfused in last
3 months?

YES

NO

Perform phenotype
for identified
antibody-antigen pair
if not already known

Refer to RCI (9)

Is patient pregnant
with an antibody
requiring quantification/
titration?

YES

NO

Perform phenotype for identified antibody-antigen pair, if not already known

It is imperative that a phenotype is performed in order to confirm that the identified antibody is an allo-antibody. E.g. patient with anti-K is K antigen negative. This will inform blood selection, should it be required. Ensure that the patient's auto result is considered prior to performing phenotyping tests. A positive auto may interfere with phenotyping tests and result in false positives. The leads to incorrectly assigning the immunological origin of an antibody as an auto-antibody which may affect blood selection. If your laboratory does not have suitable phenotyping reagents please refer to RCI for confirmation.





Is antibody screen positive?

NO

YES

Refer to RCI
(5)

Check Patient History

Does patient have a known antibody?

YES

NO

Perform Antibody ID on
Primary Panel

Has a pan-reactive antibody
been identified by IAT?

YES

NO

NO to 1 or 2

Check Patient History

1. Have antibody(ies) been identified?
2. Have additional clinically significant antibodies been excluded?
3. Is patient pregnant?
4. Has this previously been investigated by RCI?

YES to 1, 2 and 3

YES to 1 and 2

Check Results

Does patient have Rh system antibody with positive auto? Consider patient ethnicity and IBGRL referral

Refer to RCI
(18i)

Optional retrospective referral (if YES to 1&2 but NO to 4)
Refer to RCI (8)

[Click to learn more](#)

Refer to RCI (9)

If your laboratory does not have suitable phenotyping reagents, please refer to RCI.

YES

Refer to RCI (9)


Is patient pregnant with an antibody requiring quantification/titration?

YES

NO

Phenotyping results

- Patient is M negative and K negative
- Group on second sample confirmed the patient is O RhD+
- Your lab also performs extended Rh typing:





Is antibody screen positive?

NO

YES

Refer to RCI
(5)

Check Patient History

Does patient have a known
antibody?

YES

NO

Perform Antibody ID on
Primary Panel

Has a pan-reactive antibody
been identified by IAT?

YES

NO

NO to 1 or 2

Check Patient History

1. Have antibody(ies) been identified?
2. Have additional clinically significant antibodies been excluded?
3. Is patient pregnant?
4. Has this previously been investigated by RCI?

YES to 1, 2 and 3

YES to 1 and 2

Check Results

Does patient have Rh
system antibody with
positive auto? Consider
patient ethnicity and
IBGRL referral

Refer to RCI
(18i)

Optional
retrospective
referral (if YES to
1&2 but NO to 4)
Refer to RCI (8)

NO

Has patient been
transfused in last
3 months?

YES

Perform phenotype
for identified
antibody-antigen pair
if not already known

Refer to RCI (9)

Is patient pregnant
with an antibody
requiring quantification/
titration?

YES

NO





Is antibody screen positive?



YES

Refer to RCI
(5)

Check Patient History

Does patient have a known antibody?



YES

NO

Perform Antibody ID on
Primary Panel



Has a pan-reactive antibody
been identified by IAT?

YES

NO

NO to 1 or 2

Check Patient History



1. Have antibody(ies) been identified?
2. Have additional clinically significant antibodies been excluded?
3. Is patient pregnant?
4. Has this previously been investigated by RCI?

YES to 1, 2 and 3

YES to 1 and 2

Check Results



Does patient have Rh system antibody with positive auto? Consider patient ethnicity and IBGRL referral

Refer to RCI
(18i)

[Click to learn more](#)

Antibodies have been identified and all additional antibodies have been excluded

Is patient pregnant with an antibody requiring quantification/titration?


If an investigation has been successfully concluded and the patient is pregnant, consideration must be made for the need for antibody titration or quantification.


If these procedures are not performed in your laboratory, please refer to RCI.

with an antibody requiring quantification/titration?

YES

NO







Is antibody screen positive?



NO

YES

Refer to RCI
(5)

i

Check Patient History

Does patient have a known
antibody?

YES

NO

Perform Antibody ID on
Primary PanelHas a pan-reactive antibody
been identified by IAT?

YES

NO

NO to 1 or 2

Check Patient History

1. Have antibody(ies) been identified?
2. Have additional clinically significant antibodies been excluded?
3. Is patient pregnant?
4. Has this previously been investigated by RCI?



YES to 1, 2 and 3

YES to 1 and 2

Check Results

Does patient have Rh
system antibody with
positive auto? Consider
patient ethnicity and
IBGRL referralRefer to RCI
(18i)Optional
retrospective
referral (if YES to
1&2 but NO to 4)
Refer to RCI (8)Has patient been
transfused in last
3 months?

NO

YES

Perform phenotype
for identified
antibody-antigen pair
if not already known

Refer to RCI (9)

Is patient pregnant
with an antibody
requiring quantification/
titration?

YES

NO


Transfusion
Required?

NO

YES


Go to Transfusion section

Blood and Transplant



**Check Patient History**

Has patient been transfused in
the last 6 weeks?

YES**NO**

62-year-old female

Group & Screen

2 units RBC

“Anaemia”

**Not previously seen by your
lab, but transfused at
another hospital 5 years ago**

No historical antibodies

Not previously seen by RCI



Check Patient History

Has patient been transfused in
the last 6 weeks?

YES

NO





Check Patient History

Has patient been transfused in the last 6 weeks?

YES

NO



Transfusion Required?

NO

YES



62-year-old female

Group & Screen

2 units RBC


“Anaemia”

Not previously seen by your lab, but transfused at another hospital 5 years ago

No historical antibodies

Not previously seen by RCI





**Check Patient History**

Has patient been transfused in the last 6 weeks?

YES**NO**

Transfusion Required?


YES

Select red cells which are ABO compatible, D and K compatible, antigen negative for any clinically significant alloantibodies and in line with special requirements for this patient.

i**Click to learn more**

Select red cells which are **ABO compatible, D and K compatible, antigen negative for any clinically significant alloantibodies, and in line with special requirements for this patient.**

Follow your local procedures in accordance with BSH guidelines. For further information on Pre-Transfusion Compatibility Procedures in Blood Transfusion Laboratories, click [here](#).

**YES**


IAT XM compatible?

NO



BSH Guidelines – clinical significance of red cell antibodies

System	Specificity	Likely clinical significance in transfusion	Recommendation for selection of red cells for transfusion
Kell	Anti-K	Yes	Antigen negative
MNS	Anti-M (active at 37°C)	Yes	Antigen negative
Lu	Anti-Lu ^a	No	IAT crossmatch compatible at 37 °C





Check Patient History
Has patient been transfused in the last 6 weeks?

YES

NO

Transfusion Required?

NO

YES

i
Select red cells which are ABO compatible, D and K compatible, antigen negative for any clinically significant alloantibodies and in line with special requirements for this patient.

YES

IAT XM compatible?

NO

Negative reactions by IAT with patient plasma vs both O R1R1 K- M- crossmatched units





Check Patient History
Has patient been transfused in the last 6 weeks?


YES

NO

Transfusion Required?


NO


YES


 Select red cells which are ABO compatible, D and K compatible, antigen negative for any clinically significant alloantibodies and in line with special requirements for this patient.

YES

IAT XM compatible?









Case 1

- 65-year-old male
- Group & Screen
- No blood requested
- Pre-op assessment
- Historical alloanti-K
- Transfused 10 weeks ago
- LIMS and Sp-ICE results show patient as O+
E- K- with alloanti-K
- Previously seen by RCI



**Patient group &
screen results**

Case 1 – ABID Panel 1

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT	ENZ
1	R ₁ ^w R ₁	+	+	0	0	+	0	+	0	+	1	0	0	+	0	+	0	+	0	+	0		+	(+)
2	R ₁ R ₁	+	+	0	0	+	+	0	+	0	1	+	+	+	0	0	+	0	+	0	+		2	2
3	R ₂ R ₂	+	0	+	+	0	0	+	0	+	0	0	0	+	0	0	+	+	W	0	+		0	0
4	r' ^r r	0	+	0	+	+	0	+	+	0	4	0	0	+	0	0	+	0	+	+	+		0	0
5	r'' ^r r	0	0	+	+	+	+	0	+	0	3	0	0	+	0	0	+	0	+	+	0		(+)	+
6	rr	0	0	0	+	+	0	+	0	+	0	0	+	+	0	0	+	+	0	+	0	Do(a+b+)	2	2
7	rr	0	0	0	+	+	+	0	0	+	0	0	0	+	0	+	0	0	+	0	+		0	0
8	rr	0	0	0	+	+	0	+	0	+	2	0	0	+	+	+	0	0	+	+	0	HLA+	+	2
9	rr	0	0	0	+	+	0	+	+	0	3	0	0	+	0	+	0	+	0	0	+	Cob+ HLA+	0	0
10	rr	0	0	0	+	+	+	0	0	+	0	+	0	+	0	0	+	+	0	0	+		0	0
																						Auto	2	/
																						K control	2	/



Case 1 – Antibody screen

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT
1	R ₁ ^w R ₁	+	+	0	0	+	0	+	+	0	0	0	+	0	0	+	+	0	0	+		0
2	R ₂ R ₂	+	0	+	+	0	+	0	0	+	0	0	+	0	+	0	0	+	+	0		1
3	rr	0	0	0	+	+	+	0	+	0	4	+	0	0	0	+	+	0	+	0		2

DAT profile	Anti-IgG	Anti-IgA	IgM	C3c	C3d	Ctrl
Result	2	0	0	0	0	0

Case 2

- 75-year-old male
- Group & Screen
- 2 units RBC
- “Anaemia”
- Not previously seen by your lab
- Transfused at a different HTL 3 weeks ago
- Sp-ICE results show patient as A+ E- K- with no historical antibodies
- Previously seen by RCI



Patient group & screen results

Case 2 – ABID Panel 1

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT	ENZ
1	R ₁ ^w R ₁	+	+	0	0	+	0	+	0	+	0	0	0	+	0	+	0	0	+	0	+		0	0
2	R ₁ R ₁	+	+	0	0	+	+	0	+	0	4	0	+	+	0	+	0	+	0	+	0		4	2
3	R ₂ R ₂	+	0	+	+	0	0	+	0	+	3	0	0	+	0	0	+	0	+	+	0	HLA+ Do(a+b+)	2	2
4	r'r	0	+	0	+	+	0	+	0	+	0	0	0	+	0	0	+	+	0	+	+		2	2
5	r''r	0	0	+	+	+	+	0	+	0	0	0	0	+	0	0	+	+	0	0	+		2	0
6	rr	0	0	0	+	+	+	0	0	+	2	0	+	0	0	0	+	0	+	0	+		0	0
7	rr	0	0	0	+	+	0	+	0	+	1	0	+	+	0	+	0	0	+	+	0		2	2
8	rr	0	0	0	+	+	+	0	0	+	0	0	0	+	+	0	+	+	0	0	+	Cob+	0	0
9	rr	0	0	0	+	+	0	+	+	0	4	0	0	+	0	0	+	0	+	0	+		3	0
10	rr	0	0	0	+	+	+	0	0	+	4	+	0	+	0	0	+	+	0	+	0	Cob+	2	2
																						Auto	3	/
																						K control	2	/



Case 2 – Antibody screen

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	0	0	+	0	0	+	+	0	0	+		1
2	R ₂ R ₂	+	0	+	+	0	0	+	0	+	0	0	+	0	+	0	0	+	+	0		1
3	rr	0	0	0	+	+	+	0	+	0	4	+	0	0	0	+	+	0	+	0		2

DAT profile	Anti-IgG	Anti-IgA	IgM	C3c	C3d	Ctrl
Result	3	0	0	0	0	0

Case 3

- 54-year-old male
- Group & Screen
- No blood requested
- RTA
- Not previously seen by your lab
- No record on Sp-ICE
- Not previously seen by RCI
- Never been transfused



Patient group & screen results

Case 3 – ABID Panel 1

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT	ENZ
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	3	0	+	0	0	0	+	0	+	0	+		1	0
2	R ₁ R ₁	+	+	0	0	+	+	+	0	+	1	0	0	+	0	+	0	+	0	+	0	HLA+	0	0
3	R ₂ R ₂	+	0	+	+	0	0	+	0	+	3	0	0	+	0	+	0	0	+	+	0		0	0
4	r ¹ r	0	+	0	+	+	0	+	0	+	0	0	0	+	0	+	0	+	0	0	+		0	0
5	r ² r	0	0	+	+	+	+	0	+	0	1	0	0	+	0	0	+	+	0	0	+		0	0
6	rr	0	0	0	+	+	0	+	0	+	4	0	+	0	0	0	+	0	+	+	0	HLA+	0	0
7	rr	0	0	0	+	+	0	+	0	+	0	0	+	+	0	0	+	W	0	+	0		0	0
8	rr	0	0	0	+	+	+	0	0	+	0	0	0	+	+	0	+	+	0	0	+		1	0
9	rr	0	0	0	+	+	0	+	+	0	1	0	0	+	0	0	+	+	0	+	0		0	0
10	rr	0	0	0	+	+	+	0	0	+	3	+	0	+	0	0	+	0	+	+	0		0	0
																						Auto	0	/
																						K control	2	/

Case 3 – ABID Panel 1



Cell	Rh	D	C	E	c	e	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT	ENZ	Cold IAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	3	0	+	0	0	0	+	0	+	0	+		1	0	2
2	R ₁ R ₁	+	+	0	0	+	+	+	0	+	1	0	0	+	0	+	0	+	0	+	0	HLA+	0	0	1
3	R ₂ R ₂	+	0	+	+	0	0	+	0	+	3	0	0	+	0	+	0	0	+	+	0		0	0	0
4	r ¹ r	0	+	0	+	+	0	+	0	+	0	0	0	+	0	+	0	+	0	0	+		0	0	0
5	r ² r	0	0	+	+	+	+	0	+	0	1	0	0	+	0	0	+	+	0	0	+		0	0	2
6	rr	0	0	0	+	+	0	+	0	+	4	0	+	0	0	0	+	0	+	+	0	HLA+	0	0	0
7	rr	0	0	0	+	+	0	+	0	+	0	0	+	+	0	0	+	W	0	+	0		0	0	0
8	rr	0	0	0	+	+	+	0	0	+	0	0	0	+	+	0	+	+	0	0	+		1	0	2
9	rr	0	0	0	+	+	0	+	+	0	1	0	0	+	0	0	+	+	0	+	0		0	0	0
10	rr	0	0	0	+	+	+	0	0	+	3	+	0	+	0	0	+	0	+	+	0		0	0	2
																						Auto	0	/	/
																						K control	2	/	/

Case 3 – Antibody screen

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT
1	R ₁ ^w R ₁	+	+	0	0	+	0	+	+	0	0	0	+	0	0	+	+	0	0	+		0
2	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	+	0	+	0	0	+	+	0		0
3	rr	0	0	0	+	+	+	0	+	0	4	+	0	0	0	+	+	0	+	0		+/-

Case 4

- 54-year-old female
- Group & Screen
- 2 units RBC
- PV bleed
- Not previously seen by your lab
- No record on Sp-ICE
- No historical antibodies
- Not previously seen by RCI



Patient group & screen results

Case 4 – ABID Panel 1

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT	ENZ
1	R ₁ ^w R ₁	+	+	0	0	+	0	+	0	+	5	0	0	+	0	0	+	+	0	0	+		(+)	0
2	R ₁ R ₁	+	+	0	0	+	+	0	+	0	4	0	+	+	0	0	+	0	+	+	0	HLA+	0	0
3	R ₂ R ₂	+	0	+	+	0	0	+	0	+	1	0	0	+	0	+	0	+	0	0	+	HLA+	(+)	2
4	r ¹ r	0	+	0	+	+	0	+	0	+	0	0	0	+	0	+	0	+	0	0	+		(+)	0
5	r ² r	0	0	+	+	+	+	0	+	0	2	0	0	+	0	0	+	0	+	+	0		0	2
6	rr	0	0	0	+	+	+	+	0	+	3	0	+	0	0	0	0	+	0	0	+		(+)	2
7	rr	0	0	0	+	+	0	+	0	+	0	0	+	+	0	0	+	+	0	+	0	Do(a+b+)	(+)	2
8	rr	0	0	0	+	+	+	0	0	+	0	0	0	+	+	0	+	0	+	0	+		0	2
9	rr	0	0	0	+	+	0	+	+	0	2	0	0	+	+	+	0	0	+	+	0	Cob+	0	2
10	rr	0	0	0	+	+	+	0	+	0	4	+	0	+	0	0	+	+	0	0	+	Kna-	(+)	2
																						Auto	0	/
																						K control	2	/

Case 4 – Antibody screen


Cell	Rh	D	C	E	c	e	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT
1	R ₁ ^w R ₁	+	+	0	0	+	0	+	+	0	0	0	+	0	0	+	0	+	0	+		0
2	R ₂ R ₂	+	0	+	+	0	+	0	0	+	0	0	+	0	+	0	+	0	+	0		(+)
3	rr	0	0	0	+	+	+	0	+	0	4	+	0	0	0	+	+	0	0	+		+

Case 4 – secondary panel

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT	EIAT
1	R ₁ ^w R ₁	+	+	0	0	+	0	+	0	+	0	+	0	+	0	0	+	0	+	+	0		0	0
2	R ₁ R ₁	+	+	0	0	+	+	0	+	0	1	0	+	+	0	0	+	+	0	0	+		(+)	0
3	R ₂ R ₂	+	0	+	+	0	0	+	0	+	0	0	0	+	0	0	0	+	0	0	+		(+)	2
4	r' ^r r	0	+	0	+	+	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+		(+)	0
5	r'' ^r r	0	0	+	+	+	+	0	+	0	0	0	0	+	0	0	+	0	+	+	0		0	2
6	rr	0	0	0	+	+	+	+	0	+	4	0	+	0	0	0	+	+	0	+	0		(+)	2
7	rr	0	0	0	+	+	0	+	0	+	1	0	+	+	0	+	0	0	+	+	0		0	2
8	rr	0	0	0	+	+	+	0	+	0	0	0	0	+	+	0	0	0	+	0	+		0	2
9	rr	0	0	0	+	+	0	+	0	+	0	+	0	+	0	+	0	0	+	+	0		0	2
10	rr	0	0	0	+	+	0	+	+	0	3	0	0	+	0	0	0	+	0	0	+		(+)	2
																						Auto	/	/
																						K control	2	/

Case 5

- 90-year-old female
- Group & Screen
- No blood requested
- “On Dara”
- Extended phenotype results available on LIMS and Sp-ICE
- Known pan-reactive antibody due to monoclonal antibody therapy
- No alloantibodies
- DTT panel performed by RCI 2 weeks ago
- Never been transfused



Case 5 – Antibody screen

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT
1	R ₁ ^w R ₁	+	+	0	0	+	0	+	+	0	0	0	+	0	0	+	0	+	0	+		1
2	R ₂ R ₂	+	0	+	+	0	+	0	0	+	0	0	+	0	+	0	+	0	+	0		1
3	rr	0	0	0	+	+	+	0	+	0	4	+	0	0	0	+	+	0	+	0		1

Extended phenotype

D	C	E	c	e	M	N	S	s	K	k	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b
+	+	-	+	+	-	+	+	+	-	+	+	+	+	+	+	-



NHS

Blood and Transplant

